

# AI System for Detecting Flat Bladder Tumours in Cystoscopic Images

**Authors:** \*Jun Mutaguchi,<sup>1</sup> Masahiro Oda,<sup>2,3</sup> Tokiyoshi Tanegashima,<sup>1</sup> Shigehiro Tsukahara,<sup>1</sup> Shunsuke Goto,<sup>1</sup> Satoshi Kobayashi,<sup>1</sup> Takashi Matsumoto,<sup>1</sup> Masaki Shiota,<sup>1</sup> Kensaku Mori,<sup>2,3,4</sup> Masatoshi Eto<sup>1</sup>

1. Department of Urology, Kyushu University Hospital, Fukuoka, Japan
2. Information Technology Center, Nagoya University, Japan
3. Graduate School of Informatics, Nagoya University, Japan
4. Research Center for Medical Bigdata, National Institute of Informatics, Tokyo, Japan

\*Correspondence to  
mutaguchi.jun.674@m.kyushu-u.ac.jp

**Disclosure:** Olympus Corporation supported the following instruments to conduct this research (CLV-S400, OTV-S400, WM-NP2, LMD-X310S, HVO-4000ST, WA2TA430A, WA03200A, CH-S400-XZ-EB, WA05990A).

**Keywords:** AI, bladder tumour, cystoscopy, object detection, tumour detection.

**Citation:** EMJ Urol. 2025;13[1]:57-58.  
<https://doi.org/10.33590/emjurol/LVAA6547>

## INTRODUCTION AND OBJECTIVE

Bladder tumours have a high intravesical recurrence rate after transurethral resection of the bladder tumour (TURBT).<sup>1,2</sup> Overlooking tumours during TURBT occasionally causes early intravesical recurrences. In particular, flat tumours are sometimes overlooked because those tumours are more difficult to detect than elevated papillary tumours.<sup>3</sup> Detecting flat tumours is clinically beneficial because those tumours represent carcinoma *in situ* and are high-risk tumours amongst bladder tumours. Recently, AI has been developed as a new diagnosis system. You Only Look Once (YOLO) is a fast object detection AI system,<sup>4</sup> and this is a deep neural network that detects objects in an image. This technique has the potential to improve flat tumour detection rate in cystoscopy; however, only a few studies reported the utility of AI for detecting flat tumours

in cystoscopy. In this study, the authors aimed to develop an AI system to detect flat tumours in cystoscopic images.

## MATERIALS AND METHODS

The authors constructed an object detection system based on YOLO. They obtained cystoscopic images of bladder tumours during TURBT procedures at their institutions between December 2019–September 2021. Those cystoscopic images obtained were divided into training and testing data in an 8:2 ratio. YOLO was trained to detect bladder tumours in an image by learning the tumour features from training data. AI testing was conducted using YOLO, and its sensitivity, specificity, and area under the receiver operating characteristic curve were evaluated for flat tumour images.


## RESULTS

The authors obtained 2,371 tumour images and 254 non-tumour images. In tumour images, 436 flat tumour images were included. They used 1,896 images for the AI training and assessed the accuracies of the AI system. The sensitivity and specificity were 90.0% and 90.6%, respectively. With a likelihood score of 0.163, the area under the receiver operating characteristic curve in detecting flat tumours was 93.4%. An example of the detection is shown in [Figure 1](#).

## CONCLUSION

The system proposed by the authors can detect flat bladder tumours in cystoscopic images. This system has the possibility to improve detection accuracy during TURBT and might be beneficial to reduce the recurrence rate of bladder tumours after TURBT in the future.

Figure 1: Representative flat tumour detection by You Only Look Once.



The authors' proposed AI system could detect small and large flat bladder tumours. The lower right of images represents the whiteness lesion that was a carcinoma *in situ*, and their proposed AI could detect this flat lesion.

### References

1. Mutaguchi J et al. Artificial Intelligence system for detecting flat bladder tumours in cystoscopic images. Abstract A0681. EAU25, 21-24 March, 2025.
2. Kausch I et al. Photodynamic diagnosis in non-muscle-invasive bladder cancer: a systematic review and cumulative analysis of prospective studies. *Eur Urol.* 2010;57:595-606.
3. Ye Z et al. A comparison of NBI and WLI cystoscopy in detecting non-muscle-invasive bladder cancer: a prospective, randomized and multi-center study. *Sci Rep.* 2015;5:10905.
4. Redmon J, Ali F. YOLOv3: an incremental improvement. *arXiv.* 2018;DOI:1804.02767.